PrimeTime qPCR Primers

For use with PrimeTime qPCR Primers in intercalating dye assays

- 1. Centrifuge PrimeTime qPCR Primer tubes at 750 x q for 10 seconds to ensure that the primers are at the bottom of the tube.
- 2. Resuspend assay in IDTE Buffer (Cat no. 11-01-02-05) to the recommended volume, depending on the desired stock volume:

	Final desired stock concentration		
	40X	20X	10X
PrimeTime qPCR Primers, 500 rxn	250 μL	500 μL	1000 μL

Notes: During resuspension, pipette the solution up and down the sides of the tube to ensure maximum product recovery. Number of reactions is based on a $20 \, \mu L$ reaction volume.

3. Centrifuge the resuspended assay at 750 x g for 10 seconds.

Note: The resuspended primers will yield a final 1X concentration of 500 nM primers.

Storage: PrimeTime qPCR Primers are shipped dried down and can be stored up to 2 weeks at room temperature. Store resuspended assays at -20°C.

For Research Use Only. Not for use in diagnostic procedures.

© 2019 Integrated DNA Technologies, Inc., All rights reserved. PrimeTime is a trademark of Integrated DNA Technologies, Inc., and is registered in the U.S. For specific trademark and licensing information, see www.idtdna.com/trademarks. PCR-10077-PR 06/19



