Exceptional thermal stability of PrimeTime Gene Expression Master Mix enables shipping at ambient temperature, saving laboratory and environmental resources

Derek M Thomas, Maureen Young, Scott D Rose, Aurita A Menezes, Lynette A Lewis Integrated DNA Technologies, 1710 Commercial Park, Coralville, IA 52241, USA

Abstract

PrimeTime Gene Expression Master Mix is optimized for probe-based assays in two-step, reverse transcription qPCR (RT-qPCR). This master mix delivers high efficiency qPCR under fast or standard cycling conditions. Here, we show evidence of excellent benchtop stability, and demonstrate consistent qPCR performance in singleplex and duplex reactions after the master mix has been exposed to extreme temperature conditions (50°C for up to 7 days; up to 20 freeze-thaw cycles). These treatments exceed conditions experienced by the master mix during routine shipping or laboratory use. The reliable performance of the master mix allows for ambient shipping, providing environmental benefits and research dollar savings.

Introduction

Since its beginnings, IDT has striven to improve sustainability practices, and our manufacturing facilities (Coralville, IA; and San Diego, CA) have achieved ISO 14001:2004 certification for environmental management systems. Much of the focus has been on manufacturing processes to reduce air emissions, water and energy usage, and hazardous and landfill waste production. Here, we show results from experiments supporting ambient temperature shipping for PrimeTime Gene Expression Master Mix to further mitigate environmental impacts of our products.

Eliminating shipping in insulated containers with gel packs and dry ice reduces shipping costs by 35–70% (depending on your location), and may shorten delivery times. Ambient shipping also lessens detrimental environmental impacts in several ways, such as:

- Eliminating manufacturing resources required for insulated containers, dry ice, and gel packs
- Eliminating waste disposal of insulated containers or gel packs in landfills or incinerators
- Decreasing transportation resources (e.g., fuel consumption and greenhouse gas emissions) by decreasing shipping weight

See what more we can do for you at www.idtdna.com.

white paper

For consistency and convenience, many researchers use commercial master mixes that have been optimized, often with enhancers and stabilizers, to work in most standard assays. IDT scientists have developed a robust qPCR master mix for probe-based assays. The versatile PrimeTime Gene Expression Master Mix (Cat # 1055770, 1055772, or 1055771) produces consistent results under both standard and fast cycling conditions, and is compatible with a wide range of real-time PCR platforms.

We show that PrimeTime master mix is also compatible with overnight and high throughput experiments, based on the benchtop stability of the master mix. We have further extended our thermal stability testing to show that function of this master mix is not negatively affected by conditions that could be encountered during ambient shipping.

Results and discussion

High PCR efficiency after extended exposure to benchtop temperatures

We have shown that 3 batches of PrimeTime Gene Expression Master Mix produce reliable results, whether the qPCR was run immediately after setup or after the reactions were maintained at room temperature for 24 hr (Figure 1).




Figure 1. Consistent and precise results after 24 hr at room temperature using different batches of PrimeTime Gene Expression Master Mix. qPCRs consisting of the PrimeTime HPRT qPCR Assay, PrimeTime master mix, reference dye, and varying amounts of gBlocks Gene Fragments (3 of 7 dilutions shown: 10^7-10^1 copies, 8 replicates) were run immediately (Day 0) or after 24 hr (Day 1). The ΔC_q values for template levels >10 copies were <0.5. To see additional results, visit www.idtdna.com/qPCRmastermix.

We also observed similar, consistent results after testing 3 additional batches of master mix with other qPCR assays (TFRC and ACN9, B2M and IPO8, POLR2A and HMBS), run as singleplex and duplex reactions after extended benchtop time up to 72 hr (Figure 2). For all 6 assays, PCR efficiency was in the range of 84.9–103.4% ($R^2 > 0.991$) under singleplex and duplex conditions. When comparing the controls (qPCRs run immediately after reaction setup) to reactions held for 24, 48, or 72 hr at room temperature before cycling, ΔC_q values were <0.5, with a few exceptions (see Figure 2C for representative data). Together, these results indicate that there are no detectable detrimental effects from maintaining the qPCR at room temperature for up to 72 hr.

Our results show consistent qPCR results even after the reaction remained at room temperature for up to 72 hr. This highlights the stability of PrimeTime Gene Expression Master Mix and supports the use of this master mix in overnight or automated, high throughput experiments.

10²

Average $|\Delta C_q|$


0.15

0.07

-0.01

0.06

white paper



|--|

			Δ	C _q , PrimeTime	e qPCR Ass	ays		
# of copies of gBlocks fragments	HPRT	HPRT (duplexed with GUSB)	GUSB	GUSB (duplexed with HPRT)	TFRC	TFRC (duplexed with ACN9)	ACN9	ACN9 (duplexed with TFRC)
10 ⁷	-0.08	-0.21	-0.55	-0.22	-0.19	-0.13	-0.15	-0.23
104	0.04	-0.20	-0.03	-0.19	-0.07	-0.06	-0.15	-0.16
105	-0.34	-0.27	-0.22	-0.16	-0.01	-0.05	-0.15	-0.24
104	-0.10	-0.29	-0.45	-0.27	0.07	-0.02	-0.04	-0.16
10 ³	-0.05	-0.38	-0.16	-0.25	-0.08	-0.24	-0.19	-0.29
10 ²	-0.22	0.04	-0.09	-0.35	-0.40	-0.13	-0.17	0.14
Average $ \Delta C_q $	0.14	0.23	0.25	0.24	0.14	0.10	0.14	0.20
# of copies of gBlocks fragments	B2M	B2M (duplexed with IPO8)	IPO8	IPO8 (duplexed with B2M)	POLR2A	POLR2A (duplexed with HMBS)	HMBS	HMBS (duplexed with POLR2A)
107	0.01	-0.08	-0.10	-0.12	-0.07	-0.11	0.62	0.30
106	-0.01	-0.08	-0.14	-0.10	-0.06	0.04	0.07	0.15
105	-0.04	-0.06	-0.09	-0.08	-0.01	0.07	-0.01	0.14
104	-0.17	-0.09	-0.16	-0.13	0.01	0.05	0.15	0.22
10 ³	-0.01	0.03	-0.21	-0.08	-0.15	0.06	0.12	0.28

Figure 2. PrimeTime Gene Expression Master Mix has excellent benchtop stability-consistent, high PCR efficiency was observed after reactions were set up and run after 24, 48, and 72 hr at room temperature. qPCRs included PrimeTime Gene Expression Master Mix and PrimeTime qPCR Assays (ACN9, TFRC, HPRT, GUSB, B2M, IPO8, POLR2A, and HMBS) in singleplex and duplex. The qPCRs were set up and either run immediately (control) or remained at room temperature for 24, 48, or 72 hr before cycling. Consistent, high efficiency PCR results are shown by representative amplification curves (A) and standard curves (B) from ACN9 assays run in singleplex: controls (gray) and 24, 48, and 72 hr samples (blue, green, and orange, respectively). (C) Representative data comparing ΔC_q values of the control qPCRs to reactions held for 72 hr at room temperature show that the difference in the average ΔC_q values from technical triplicates were within ±0.50, except in 3 instances (orange).

-0.28

0.13

0.22

-0.13

0.07

0.03

0.06

0.01

0.16

0.34


0.24

Consistent C_{α} values after extended heat-stress at 50°C

To further characterize the performance dependability and consistency of PrimeTime master mix, we examined data from accelerated stability experiments. The extreme conditions, which would not be encountered during normal product life and use, included elevated temperatures for extended periods of time.

Three batches of master mix were heated at 50°C for 1, 3, or 7 days before use in qPCR experiments (Figure 3). We chose HPRT and GUSB as our representative qPCR assays for this analysis. High PCR efficiency was observed (range of 95.1–103.6% with R² >0.991) with all assays tested. When comparing the controls (containing frozen, unheated master mix) to reactions containing heat-treated master mix, ΔC_q values were <0.5, with a few exceptions that were primarily seen when using a low amount of starting template (see Figure 3C for representative data). The small ΔC_q values emphasize the stability and robustness of PrimeTime master mix after extreme thermal stress treatment.

Our results consistently show successful qPCR results even after PrimeTime Gene Expression Master Mix was heated at 50°C for 7 days, which emphasizes the stability of the master mix and its components, when not stored at -20°C.



$C \Delta C_q = Aver$	rage C _q (con	trols) – Avera	ige C _q (X day	s)				
-	∆C _q , HPRT			ΔC_q , HPRT (duplexed with GUSB)				
# of copies of gBlocks fragments	Control vs. 1 day at 50° C	Control vs. 3 days at 50° C	Control vs. 5 days at 50°C	Control vs. 7 days at 50°C	Control vs. 1 day at 50°C	Control vs. 3 days at 50°C	Control vs. 5 days at 50°C	Control vs. 7 days at 50°C
107	-0.06	-0.15	-0.02	-0.01	-0.11	-0.09	0.06	0.03
106	-0.07	-0.14	-0.06	0.03	-0.10	-0.13	-0.04	0.12
10 ⁵	-0.05	-0.16	-0.02	0.05	-0.09	-0.16	0.03	0.05
104	-0.08	-0.14	-0.03	0.06	-0.06	-0.02	0.14	0.11
10 ³	0.05	-0.12	-0.01	0.03	-0.01	0.00	0.01	0.14
10 ²	-0.07	-0.08	-0.04	-0.03	-0.11	-0.02	0.00	0.00
Average ΔC _q	0.06	0.13	0.03	0.04	0.08	0.07	0.05	0.08
	ΔC _q , GUSB				ΔC _q , GUSB (duplexed with HPRT)			
# of copies of gBlocks fragments	Control vs. 1 day at 50° C	Control vs. 3 days at 50° C	Control vs. 5 days at 50°C	Control vs. 7 days at 50°C	Control vs. 1 day at 50°C	Control vs. 3 days at 50°C	Control vs. 5 days at 50°C	Control vs. 7 days at 50°C
107	-0.06	-0.04	0.04	0.10	-0.06	-0.06	0.13	0.01
106	0.20	0.19	0.37	0.35	0.08	0.03	0.21	0.27
10 ⁵	-0.03	0.03	0.12	0.09	-0.02	-0.03	0.15	0.08
104	-0.02	-0.02	0.18	0.14	0.05	0.05	0.22	0.20
10 ³	0.13	0.18	0.27	0.25	0.04	-0.13	0.13	0.03
10 ²	0.10	0.28	0.08	0.29	0.18	0.13	0.22	0.06
Average $ \Delta C_q $	0.09	0.12	0.18	0.20	0.07	0.07	0.18	0.11

Figure 3. Even after elevated temperature treatments (50°C for up to 7 days), PrimeTime Gene Expression Master Mix provides high PCR efficiency with no change in C_q values, when compared to untreated master mix. PrimeTime master mix was incubated at 50°C for 1, 3, 5, or 7 days (experimental), or stored at –20°C until use (control). qPCRs included the heat-treated or control master mix and PrimeTime qPCR Assays (HPRT and GUSB) in singleplex and duplex. Consistent, high efficiency PCR results are shown by representative amplification curves (A) and standard curves (B) from HPRT assays (singleplex) that remained at room temperature for 24 hr before running the thermal cycler: control master mix (gray) and heat-treated master mix (1 day, light blue; 3 days, dark blue; 5 days, green; and 7 days, orange). (C) Representative data comparin ΔC_q values of the control qPCRs to reactions that included 1-, 3-, 5-, or 7-day, heat-treated master mix show that the difference in the average C_q values from technical triplicates were within ±0.50.


Consistent C_q values after multiple freeze-thaw cycles

We also tested the effect of repeated freezing and thawing on the functional performance of PrimeTime master mix. Freeze-thawing can occur during normal use in the laboratory or, potentially, during shipping. We tested up to 20 freeze-thaw cycles, which is more extreme than typical laboratory usage.

PrimeTime Gene Expression Master Mix was frozen and thawed 5, 10, or 20 times before use in qPCR experiments (Figure 4). We chose HPRT and GUSB as our representative qPCR assays for this analysis. High PCR efficiency was observed (range of 95.4–104.9% with R² >0.994) with all qPCRs tested. When comparing the control qPCRs (containing frozen, untreated master mix) to reactions containing frozen-thawed master mix, ΔC_q values were all <0.5, with a few exceptions (see Figure 4C for representative data). The small ΔC_q values indicate essentially no detrimental effects from repeated freeze-thawing of PrimeTime master mix.

white paper

These results show that efficiency of the qPCR was unaffected, even after the PrimeTime master mix was subjected to as many as 20 freeze-thaw cycles, which further emphasizes the stability of the master mix and its components.



 ΔC_q = Average C_q (controls) – Average C_q (X freeze-thaw cycles)

	HPRT		HPRT (duplexed with GUSB)			
Control vs. 5 freeze-thaw cycles	Control vs. 10 freeze-thaw cycles	Control vs. 20 freeze-thaw cycles	Control vs. 5 freeze-thaw cycles	Control vs. 10 freeze-thaw cycles	Control vs. 20 freeze-thaw cycles	
0.49	0.13	0.08	-0.01	-0.04	0.04	
0.56	0.13	0.15	0.12	-0.10	-0.13	
0.01	0.03	-0.06	-0.34	-0.10	-0.13	
0.32	-0.18	0.14	0.21	-0.19	-0.09	
0.47	0.07	0.09	0.44	-0.21	0.28	
0.53	-0.33	0.32	-0.12	-0.22	-0.26	
0.40	0.14	0.14	0.21	0.14	0.15	
	GUSB		GUSB	(duplexed with	HPRT)	
	5 freeze-thaw cycles 0.49 0.56 0.01 0.32 0.47 0.53	Control vs. Control vs. 5 freeze-thaw 10 freeze-thaw cycles cycles 0.49 0.13 0.56 0.13 0.01 0.03 0.32 -0.18 0.47 0.07 0.53 -0.33	Control vs. Control vs. Control vs. 5 freeze-thaw 10 freeze-thaw 20 freeze-thaw 0.49 0.13 0.08 0.56 0.13 0.15 0.01 0.03 -0.06 0.32 -0.18 0.14 0.47 0.07 0.09 0.53 -0.33 0.32	Control vs. Stresset and transformed and transform	Control vs.	

# of copies of gBlocks fragments	Control vs. 5 freeze-thaw cycles	Control vs. 10 freeze-thaw cycles	Control vs. 20 freeze-thaw cycles	Control vs. 5 freeze-thaw cycles	Control vs. 10 freeze-thaw cycles	Control vs. 20 freeze-thaw cycles
107	0.12	-0.08	0.09	0.02	-0.07	0.05
106	0.13	0.13	0.14	0.07	-0.05	-0.04
105	0.04	0.12	0.08	-0.31	-0.06	-0.02
104	0.46	0.02	0.23	0.22	-0.19	-0.02
10 ³	0.16	0.11	0.27	0.44	-0.19	0.31
10 ²	0.19	0.05	0.00	0.18	-0.06	0.02
Average $ \Delta C_q $	0.18	0.09	0.14	0.21	0.11	0.08

Figure 4. Multiple freeze-thaw cycles (up to 20) does not affect the stability of the PrimeTime Gene Expression Master Mix, which still performs with consistent C_q values and high PCR efficiency. PrimeTime master mix underwent 5, 10, or 20 freeze-thaw cycles, while control master mix was stored frozen (-20°C) until use. qPCRs included the treated or control master mix and PrimeTime qPCR Assays (HPRT and GUSB) in singleplex and duplex. Consistent, high efficiency PCR results are shown by representative amplification curves (**A**) and standard curves (**B**) from HPRT assays (singleplex) that remained at room temperature for 24 hr before running the thermal cycler: control master mix (gray) and the treated master mix (5 freeze-thaw cycles, blue; 10 cycles, green; 20 cycles, orange). (**C**) Representative data comparing ΔC_q values of the control qPCRs to reactions that included frozen and thawed master mix show that the difference in the average C_q values from technical triplicates were within ±0.50, except in 2 instances (orange). F-T = freeze-thaw.

white paper

Case study

Thus far, the results from our experiments indicate that all the components of the master mix, including the hot-start polymerase, are stable under varying temperature conditions. We have tested 2 batches of PrimeTime Gene Expression Master Mix that were either stored at –20°C (control) or shipped at ambient temperatures. Roundtrip shipments from IDT facilities in the USA to Singapore or Australia took 7 or 8 days, respectively. Data about ambient conditions were available for the first 2 days of transit to Singapore. These records indicated that the temperature varied between 14 and 24°C and pressure varied between 800 and 1050 mBar.

We chose HPRT and GUSB as representative qPCR assays to compare control master mix (held at a constant –20°C) and experimental master mix (shipped at ambient temperature). High PCR efficiency was observed (range of 97–101%, R² >0.99) with all assays tested. When comparing the control reactions containing frozen, unshipped master mix to reactions containing shipped master mix, ΔC_q values were all <0.50. The small ΔC_q values indicate that the performance of the master mix was not compromised by ambient shipping.

Environmental benefits of shipping at ambient temperature vs. using dry ice

To illustrate the environmental impact of switching to ambient shipping, consider that each package shipped with dry ice carries an extra 1.5–6 kg (3–13 lbs) of gross weight. In a typical year for IDT, this can mean a decrease of ~30,400 kg (~67,000 lbs) in total package weight if 10,000 orders were shipped at ambient temperature, instead of on dry ice. Also, assuming cardboard packaging is recycled, ambient shipping materials would not create landfill waste, while dry-ice shipping materials include 28.5 tons of dry ice and would create ~437 m² (~4700 ft²) of landfill waste from the expanded polystyrene (foam) containers.

Conclusions

PrimeTime Gene Expression Master Mix has proven to be stable under various temperature conditions, which supports shipping at ambient temperature, as well as use in overnight or high throughput experiments. Consistently high PCR efficiencies and reproducible C_q results are observed even under the following conditions:

- The qPCR, including the master mix, was set up and remained at room temperatures for up to 72 hr.
- The master mix was stored at elevated temperatures for extended periods of time (up to 7 days at 50°C).
- The master mix was subjected to 20 freeze-thaw cycles.

Using the versatile PrimeTime master mix in your gene expression studies allows you to focus on your results instead of troubleshooting your qPCR conditions.

Ambient shipping saves packaging and transportation resources, and greatly reduces solid waste. Because our experiments show that PrimeTime Gene Expression Master Mix is not adversely affected by shipping at ambient conditions, we have moved to shipping at ambient temperature, which is better for your budget and the environment.

Visit www.idtdna.com/qPCRMasterMix for more information, or contact applicationsupport@idtdna.com with questions. To learn more about IDT sustainability initiatives, visit www.idtdna.com/sustainability.

Methods

Heat stressing the master mix

Vials of PrimeTime Gene Expression Master Mix from 3 batches were heated in a 50°C incubator for 1, 3, and 7 days. One of the 3 batches was incubated at 50°C for 14 days with similar results (data not shown). The unheated control master mix was stored at –20°C.

Freeze-thawing the master mix

PrimeTime Gene Expression Master Mix (1.5 mL vials) were subjected to 5, 10, and 20 freeze-thaw cycles. Each cycle consisted of thawing the master mix at room temperature for at least 1 hr, followed by freezing at –20°C for at least 12 hr. The unheated control master mix was stored at –20°C.

qPCR

All experiments shown were run on an Applied Biosystems® 7900HT Real-Time PCR System (Thermo Fisher Scientific) under fast cycling conditions (Table 1). Eight replicate qPCRs run as singleplex or duplexed assays were used in the experiments for Figure 1, while triplicate qPCRs were run as singleplex or duplexed assays for the remaining experiments. Details of the qPCRs are provided in Tables 1 and 2. For some experiments, the singleplex or duplex PCR mix was left on the benchtop (room temperature, 15–25°C) for 24, 48, or 72 hr before running the assays on the thermal cycler (Figure 2). For other experiments, the master mix was subjected to temperature stress (50°C heating or multiple freeze-thaw cycles) before use, where the PCRs were run immediately after set up (data not shown) or after 24 hr at room temperature (Figures 3 and 4).

gBlocks Gene Fragments, which are sequence verified, double-stranded DNA, were used as the qPCR template. Ten-fold dilutions were made from 10⁷ to 10 copies per reaction. Some experiments were also performed with cDNA prepared from a commercially available Universal RNA using the AffinityScript[®] RT-PCR Kit (Agilent Technologies). The results with cDNA template (5-fold dilutions from 50 to 0.0032 ng) were similar to those with gBlocks Fragment templates (data not shown).

 ΔC_q was calculated from the average C_qs of replicate reactions for the conditions being studied (e.g., qPCR results from different time points or from different temperature storage conditions).

Table 1. qPCR setup. PrimeTime Gene Expression Master Mix includes a separate tube of reference dye for use as needed. IDT offers Nuclease-Free Water (Cat # 11-04-02-01). Fast cycling conditions: 3 min. 95°C; 49 x (5 sec. 95°C, 30 sec. 60°C).

	Singleplex	Singleplex	Duplex
Assay 1—predesigned PrimeTime primers and FAM-labeled probe (1X final)	0.3 µL		0.3 µL
Assay 2—predesigned PrimeTime primers and HEX-labeled probe (1X final)	_	0.3 µL	0.3 µL
PrimeTime Gene Expression Master Mix (1X final)	6 µL	6 µL	6 µL
Reference dye (0.5 µM final)	0.24 µL	0.24 µL	0.24 µL
Nuclease-Free Water	0.46 µL	0.46 µL	0.16 µL
Template (gBlocks Gene Fragments) (101–107 copies)	5 µL	5 µL	5 µL
Reaction volume	12 µL	12 µL	12 µL

white paper

Table 2. qPCR assay and template details.

Component	Detailed information							
Primers and Probes	PrimeTime qPCR Assays: visit www.idtdna.com/PrimeTime for more information. The following pairs of predesigned assays were run in singleplex or duplex:							
	Target	Probe flourophore and quenchers	Assay ID					
	HPRT	FAM/ZEN/IBFQ	Hs.PT.58v.45621572					
	GUSB	HEX/ZEN/IBFQ	Hs.PT.58v.27737538					
	TFRC	FAM/ZEN/IBFQ	Hs.PT.58.3164874					
	ACN9	HEX/ZEN/IBFQ	Hs.PT.58.588601					
	B2M	FAM/ZEN/IBFQ	Hs.PT.58v.18759587					
	IPO8	HEX/ZEN/IBFQ	Hs.PT.58.3371286					
	POLR2A	FAM/ZEN/IBFQ	Hs.PT.58.25515089					
	HMBS	HEX/ZEN/IBFQ	Hs.PT.58.40437381					
Templates	gBlocks Gene Fragments: Visi	t www.idtdna.com/gBlocks for more inform	ation.					
	HPRT, 200 bp (5′→3′):							
	TGTTGGATTT GAAATTCCAG	ACAAGTTTGT TGTAGGATAT GCCCTTGACT	ATAATGAATA CTTCAGGGAT TTGAATCAT					
	TTTGTGTCAT TAGTGAAACT TCCTATTGAC ATCGCTTGAC	GGAAAAGCAA AATACAAAGC CTAAGATGAG AAATTATCAA TGTTCTAGTT	AGTTCAAGTT GAGTTTGGAA ACATCTGGA					
	GUSB, 170 bp (5′→3′):							
	ACTGGTATAA GAAGTATCAG	AAGCCCATTA TTCAGAGCGA GTATGGAGCA	GAAACGATTG CAGGGTTTCA CCAGGATCC					
	CCTCTGATGT TCACTGAAGA CGTGGTTGGA	GTACCAGAAA AGTCTGCTAG AGCAGTACCA	TCTGGGTCTG GATCAAAAAC GCAGAAAAT					
	TFRC, 160 bp (5'→3'):							
		CATTTGTGAG GGATCTGAAC CAATACAGAG GGAGACTTCT TCCGTGCTAC TTCCAGACTA						
	ACN9, 160 bp (5'→3'):							
		ACGAATTTAG GAGACATAAG ACCGTTGGTT TTATTGCAAC AGGCTAACGA AAACAGACAA						
	B2M, 170 bp (5'→3'):							
		TCTTTCTATC TCTTGTACTA CACTGAATTC GTCACAGCCC AAGATAGTTA AGTGGGATCG						
	IPO8, 200 bp (5'→3'):							

ATGTGGCAGC TTCTAGGTAT ACTATATGAA GTGTTTCAGC AGGATTGCTT TGAATACTTT ACAGACATGA TGCCTCTCCT GCATAATTAT GTGACAATAG ATACAGATAC CTTACTATCA AATGCAAAAC ATTTAGAAAT TCTTTTACA ATGTGTAGGA AGGTACTATG TGGAGATGCA GGAGAAGATG CAGAGTGTCA

POLR2A, 160 bp (5'→3'):

GATGACAATG CAGAGAAGCT GGTGCTCCGT ATTCGCATCA TGAACAGCGA TGAGAACAAG ATGCAAGAGG AGGAAGAGGT GGTGGACAAG ATGGATGATG ATGTCTTCCT GCGCTGCATC GAGTCCAACA TGCTGACAGA TATGACCCTG CAGGGCATCG

HMBS, 180 bp (5'→3'):

ACTCCTTGAA GGACCTGCCC ACTGTGCTTC CTCCTGGCTT CACCATCGGA GCCATCTGCA AGCGGGAAAA CCCTCATGAT GCTGTTGTCT TTCACCCAAA ATTTGTTGGG AAGACCCTAG AAACCCTGCC AGAGAAGAGT GTGGTGGGAA GTGGTGGGAA GCGAAGAGCA GCCCAGCTGC

Revision history

Version	Date released	Description of changes
1	May 2016	Original white paper
1.1	February 2017	Updated with new IDT logo and styles
2	April 2019	Corrected sequence of the HPRT template on page 10

Integrated DNA Technologies, Inc. (IDT) is your Advocate for the Genomics Age. For more than 30 years, IDT's innovative tools and solutions for genomics applications have been driving advances that inspire scientists to dream big and achieve their next breakthroughs. IDT develops, manufactures, and markets nucleic acid products that support the life sciences industry in the areas of academic and commercial research, agriculture, medical diagnostics, and pharmaceutical development. We have a global reach with personalized customer service. See what more we can do for **you** at **www.idtdna.com**.

Technical support: applicationsupport@idtdna.com

For Research Use Only. Not for use in diagnostic procedures.

© 2019 Integrated DNA Technologies, Inc. All rights reserved. gBlocks and PrimeTime are trademarks of Integrated DNA Technologies, Inc., and are registered in the USA. All other marks are the property of their respective owners. For specific trademark and licensing information, see www.idtdna.com/trademarks. PCR-10032-PR 04/2019